

Technical Name: DENTAL IMPLANTS (OSSEOINTEGRATED)

Commercial Name: IMPLANT

Trading Model:

Cylindrical Biodent HEX (all the models)

Tapered Bioneck TRI (al the models)

Cylindrical Dynamic CMI (all the models)

Tapered Single CMI (all the models)

Kort Hex (all the models)

Kort CMI Plus (todos os modelos)

Cylindrical Exakort HI (todos os modelos)

PRODUCT DESCRIPTION:

The **DÉRIG Dental Implants** are constituted: of a cylindrical or tapered body, with rough surface. This chemical treatment process increases its surface of contact, thus increasing the primary stability and accelerating the osseointegration. The Implant can have a milled area in the apex, distributed in a seamless manner in order to facilitate its insertion. In the longitudinal shaft of the upper part of the implant there is a central thread that allows fixing of the prosthetic component using a screw. It is produced as from turning of **4 degree pure titanium bars in compliance with ASTM F67**.

Os DÉRIG Implants are supplied with Cover Screw, to protect the inner parts during the osseointegration, if the surgeon decides to install a temporary prosthesis after this process.

In order to provide an option to the surgeon we have available 7 groups of implants that differ between them in the following aspects:

- a) **External profile:** cylindrical or tapered.
- b) **Thread pitch:** 0.5 to 1.2mm.
- c) **Thread profile:** Thread in V, Trapezoidal Thread and Square Thread.
- d) **External profile size:** Ø 3.0 to 5.0mm and length of 5.5 through 18mm.

- e) **Internal thread:** Ø 1.6, 1.8, 2.0 and 2.5mm.
- f) **Prosthetic connective platform:** External hexagon, internal triangle and internal tapered.
- g) **Codes:** The NP, RP, WP and UP codes refer to the size of the implant platform

Connection:

HEX – External Hexagon

TRI – Internal Triple Channel

CMI – Cone Morse Indexed (indexed morse cone)

HI – Internal Hexagon

Platform:

NP – Narrow Platform → receives NP components

RP – Regular Platform → receives RP-UP components

UP – Unified Platform → receives RP-UP components


WP – Wide Platform → receives WP components

Platform Code	Measures of platform diameters according to the system(mm)						
	Biodent HEX	Bioneck TRI	Dynami CMI	Singular CMI	Kort HEX	Kort CMI Plus	Exakort HI
NP	Ø3.5	Ø3.5	Ø2.75	Ø2.75	Ø3.5	Ø2,75	Ø3,40
RP	Ø4.1	Ø4.3	Ø2.75	Ø2.75	Ø4.1	-	-
UP	Ø4.1	Ø4.3	-	-	-	-	-
WP	Ø5.1	Ø5.0	-	-	-	-	-

Example:

Bioneck Tapered Implant TRI NP 4.3 x 10mm

Biodent Cylindrical Implant HEX RP 3.75 x 13mm



The DÉRIG Implants are divided below in 7 groups, and each group recommended for one or more regions of the dental arch due to the technique or surgical situation

adequate to the surgical plan of each patient, in order to provide options to the surgeon.

- 1) **Cylindrical Biodent HEX:** Cylindrical external profile, with double thread in “V”, with chambers in its apex, internal with Ø 1.6, 2.0 and 2.5mm thread, external prosthetic connective platform in a hexagonal format.
- 2) **Tapered Bioneck TRI:** Tapered external profile with trapezoidal thread, internal with Ø 1.8 and 2.0mm thread, prosthetic connective platform in triangular format.
- 3) **Cylindrical Dynamic CMI:** Cylindrical external Profile, with square double thread, internal with Ø 1.6mm thread, internal prosthetic connective platform, in tapered format.
- 4) **Tapered Singular CMI:** Tapered external profile with trapezoidal thread, internal with Ø 1.6mm, thread, internal prosthetic connective platform in tapered format.
- 5) **Kort HEX:** Cylindrical external profile, with double thread in “V”, with chambers in its apex and in the coronary portion, internal with Ø 1.6, and 2.0mm thread, external prosthetic connective platform in a hexagonal format.
- 6) **Kort CMI Plus:** External profile cylindrical, with double "V" thread, with chambers at its apex and in the coronary portion, internal with Ø 1,6 thread, internal prosthetic connective platform, in tapered format.
- 7) **Cylindrical Exakort HI:** Cylindrical external Profile, with square double thread, internal with M 1.8 (Ø 1.8mm) thread, internal prosthetic connective platform, in hexagonal format with tapered contact.

CYLINDRICAL BIODENT HEX:

- **Ø 3,3mm NP**, available in the following lengths: 8.5 – 10.0 – 11.5 – 13.0 – 15.0 – 16.0 and 18.0mm, Ø 1.6mm internal thread, with self-threading chambers in its apex. Connective Platform with HEX NP prosthetic components.
- **Ø 3.75mm NP**, available in the following lengths: 8.5 – 10.0 – 11.5 – 13.0 – 15.0 – 16.0 and 18.0mm, Ø 2.0mm internal thread, with self-threading

chambers in its apex. Connective Platform with HEX RP-UP prosthetic components.

- **Ø 4.0mm RP**, available in the following lengths: 8.5 – 10.0 – 11.5 – 13.0 – 15.0 – 16.0 and 18.0mm, Ø 2.0mm internal thread, with self-threading chambers in its apex. Connective Platform with HEX RP-UP prosthetic components.
- **Ø 5.0mm UP**, available in the following lengths: 8.5 – 10.0 – 11.5 – 13.0 – 15.0 – 16.0 and 18.0mm, Ø 2.0mm internal thread, with self-threading chambers in its apex. Connective Platform with HEX RP-UP prosthetic components.
- **Ø 5.0mm WP**, available in the following lengths: 8.5 – 10.0 – 11.5 – 13.0 – 15.0 – 16.0 and 18.0mm, Ø 2.5mm internal thread, with self-threading chambers in its apex. Connective Platform with HEX WP prosthetic components.

TAPERED BIONECK TRI:

- **Ø 3.5mm NP**, available in the following lengths: 8.0 – 10.0 – 11.5 – 13.0 – 15.0 – 16.0mm, Ø 1.8mm internal thread, with self-threading chambers in its apex. Connective Platform with TRI NP prosthetic components.
- **Ø 4.3mm NP**, available in the following lengths: 8.0 – 10.0 – 11.5 – 13.0 – 15.0 – 16.0mm, Ø 1.8mm internal thread. Connective Platform with DERIG TRI NP prosthetic components.
- **Ø 4.3mm RP**, available in the following lengths: 8.0 – 10.0 – 11.5 – 13.0 – 15.0 – 16.0mm, Ø 2.0mm internal thread, with self-threading chambers in its apex. Connective Platform with TRI RP-UP prosthetic components.
- **Ø 5.0 mm UP**, available in the following lengths: 8.0 – 10.0 – 11.5 – 13.0 – 15.0 – 16.0mm, Ø 2.0mm internal thread, with self-threading chambers in its apex. Connective Platform with TRI RP-UP prosthetic components.
- **Ø 5.0 mm WP**, available in the following lengths: 8.0 – 10.0 – 11.5 – 13.0 – 15.0 – 16.0mm, Ø 2.0mm internal thread, with self-threading chambers in its apex. Connective Platform with TRI WP prosthetic components.

CYLINDRICAL DYNAMIC CMI:

- **Ø 3.0mm**, available in the following lengths: 10.0 – 11.5 – 13.0 – 15.0mm, Ø 1.6mm internal thread, with self-threading chambers in its apex. Connective Platform with CMI 3.0 prosthetic components.
- **Ø 3.5mm NP**, available in the following lengths: 8.5 – 10.0 – 11.5 – 13.0 – 15.0 – 16.0, 18.0mm, Ø 1.6mm internal thread, with self-threading chambers in its apex. Connective Platform with CMI NP prosthetic components.
- **Ø 4.3mm RP**, available in the following lengths: 8.5 – 10.0 – 11.5 – 13.0 – 15.0 – 16.0, 18.0mm, Ø 1.6mm internal thread, with self-threading chambers in its apex. Connective Platform with CMI RP prosthetic components.
- **Ø 5.0 mm RP**, available in the following lengths: 8.5 – 10.0 – 11.5 – 13.0 – 15.0 – 16.0, 18.0mm, Ø 1.6mm internal thread, with self-threading chambers in its apex. Connective Platform with CMI RP prosthetic components.

TAPERED SINGLE CMI:

- **Ø 3.5mm NP**, available in the following lengths: 8.0 – 10.0 – 11.5 – 13.0 – 15.0 – 16.0mm, Ø 1,6mm internal thread. Connective Platform with CMI NP prosthetic components.
- **Ø 4.3mm RP**, available in the following lengths: 8.0 – 10.0 – 11.5 – 13.0 – 15.0 – 16.0mm, Ø 1,6mm internal thread. Connective Platform with CMI RP prosthetic components.
- **Ø 5.0 mm RP**, available in the following lengths: 8.0 – 10.0 – 11.5 – 13.0 – 15.0 – 16.0mm, Ø 1,6mm internal thread. Connective Platform with CMI RP prosthetic components.

KORT HEX:

- **Ø 4.0mm NP**, available in the following lengths: 5.5 –7.0mm, Ø 1.6mm internal thread, with self-threading chambers in its apex. Connective Platform with HEX NP prosthetic components.

- **Ø 5.0mm RP**, available in the following lengths: 5.5 - 7.0mm, Ø 2.0mm internal thread, with self-threading chambers in its apex. Connective Platform with HEX RP prosthetic components.

KORT CMI Plus:





- **Ø 3,5mm NP**, available in lengths: 5.5 - 7.0mm, internal thread Ø 1.6mm, with self-tapping chambers at its apex. Connective platform with prosthetic components CMI NP.
- **Ø 4.3mm NP**, available in lengths: 5.5 - 7.0mm, internal thread Ø 1.6mm, with self-threading chambers at its apex. Connective platform with prosthetic components CMI NP.
- **Ø 5,0mm NP**, available in lengths: 5.5 - 7.0mm, internal thread Ø 1.6mm, with self-tapping chambers at its apex. Connective platform with prosthetic components CMI NP.

CYLINDRICAL EXAKORT HI:

- **Ø 3.3mm NP**, available in the following lengths: 8.0 – 10.0 – 11.5 – 13.0 – 15.0, M1.8 (Ø 1.8mm) internal thread, with self-threading chambers in its apex. Connective Platform with HI NP prosthetic components.
- **Ø 3.8mm NP**, available in the following lengths: 8.0 – 10.0 – 11.5 – 13.0 – 15.0, M1.8 (Ø 1.8mm) internal thread, with self-threading chambers in its apex. Connective Platform with HI NP prosthetic components.
- **Ø 4.2mm NP**, available in the following lengths: 8.0 – 10.0 – 11.5 – 13.0 – 15.0, M1.8 (Ø 1.8mm) internal thread, with self-threading chambers in its apex. Connective Platform with HI NP prosthetic components.





CYLINDRICAL BIODENT HEX




Connective Platform	Platform diameter	Diameter of the implant	Lengths
	Ø3, 5 mm	Ø3, 3 mm	8, 5 mm, 10 mm, 11, 5 mm, 13 mm, 15 mm, 16 mm and 18 mm
	Ø4, 1 mm	Ø3, 75 mm	8, 5 mm, 10 mm, 11, 5 mm, 13 mm, 15 mm, 16 mm and 18 mm
		Ø4 mm	8, 5 mm, 10 mm, 11, 5 mm, 13 mm, 15 mm, 16 mm, 17 mm and 18 mm
	Ø4, 1 mm	Ø5, 0 mm	8, 5 mm, 10 mm, 11, 5 mm, 13 mm, 15 mm, 16 mm and 18 mm
	Ø5, 1 mm	Ø5, 0 mm	8, 5 mm, 10 mm, 11, 5 mm, 13 mm, 15 mm, 16 mm and 18 mm




TAPERED BIONECK TRI



Connective Platform	Platform diameter	Diameter of the implant	Lengths
	Ø3, 5 mm	Ø3, 75 mm	8 mm, 10 mm, 11, 5 mm, 13 mm, 15 mm and 16 mm
		Ø4, 3 mm	8 mm, 10 mm, 11, 5 mm, 13 mm, 15 mm and 16 mm
	Ø4, 3 mm	Ø4, 3 mm	8 mm, 10 mm, 11, 5 mm, 13 mm, 15 mm and 16 mm
	Ø4, 3 mm	Ø5, 0 mm	8 mm, 10 mm, 11, 5 mm, 13 mm, 15 mm and 16 mm
	Ø5, 0 mm	Ø5, 0 mm	8 mm, 10 mm, 11, 5 mm, 13 mm, 15 mm and 16 mm




DYNAMIC CYLINDRICAL CMI




Connective Platform	Platform diameter	Diameter of the implant	Interface with the pillar	Lengths
	Ø3, 0 mm	Ø3, 6 mm	Ø2, 75 mm	8, 5 mm, 10 mm, 11, 5 mm, 13 mm and 15 mm
	Ø3, 0 mm	Ø4, 3 mm	Ø2, 75 mm	8, 5 mm, 10 mm, 11, 5 mm, 13 mm and 15 mm
	Ø3, 0 mm	Ø5, 0 mm	Ø2, 75 mm	8, 5 mm, 10 mm, 11, 5 mm, 13 mm and 15 mm



TAPERED SINGLE CMI




Platform Connective	Platform diameter	Diameter of the implant	Interface with the pillar	Lengths
	Ø3, 0 mm	Ø3, 75 mm	Ø2, 75 mm	8 mm, 10 mm, 11, 5 mm, 13 mm and 16 mm
	Ø3, 0 mm	Ø4, 3 mm	Ø2, 75 mm	8 mm, 10 mm, 11, 5 mm, 13 mm and 16 mm
	Ø3, 0 mm	Ø5, 0 mm	Ø2, 75 mm	8 mm, 10 mm, 11, 5 mm, 13 mm and 16 mm


KORT HEX



Platform Connective	Platform diameter	Diameter of the implant	Lengths
	Ø3, 5 mm	Ø4 mm	5, 5 mm and 7 mm
	Ø4, 1 mm	Ø5, 0 mm	5, 5 mm and 7 mm


KORT CMI PLUS




Platform Connective	Interface with the pillar	Diameter of the implant	Lengths
	Ø2, 75 mm	Ø3, 5 mm	5, 5 mm and 7 mm
	Ø2, 75 mm	Ø4, 3 mm	5, 5 mm and 7 mm
	Ø2, 75 mm	Ø5, 0 mm	5, 5 mm and 7 mm

HI EXAKORT IMPLANT




Platform Connective	Interface with the pillar	Diameter of the implant	Lengths
	Ø3, 4 mm	Ø3, 3 mm	8 mm, 10 mm, 11, 5 mm, 13 mm and 15 mm
	Ø3, 7 mm	Ø3, 8 mm	8 mm, 10 mm, 11, 5 mm, 13 mm and 15 mm
	Ø3, 7 mm	Ø 4, 2 mm	8 mm, 10 mm, 11, 5 mm, 13 mm and 15 mm


All Dérig implants have internal thread with "V" profile and can vary between $\varnothing 1,60 / \varnothing 1,80 / \varnothing 2,0 / \varnothing 2,5$




Implant body with conical profile
(Bioneck and singular implants)




Biodent Implants have external thread with triangular profile



Bioneck and Singular Implants have external thread with trapezoidal profile




Dynamic implants have external thread with square profile



The implants of the Biodent line have external prosthetic connective platform in hexagonal format



The implants of the Biomeck line have external prosthetic connective platform in triangular format



The implants of the CMI line have external prosthetic connective platform in conical format



LIST OF ACCESSORIES DESIGNED TO INTEGRATE THE PRODUCT:

The **DÉRIG** Implants have the accessories described below:

COVER SCREW:

Each implant has its respective Cover Screw, also produced in degree 4 pure titanium, in compliance with ASTM F67, to be installed on the Implant, with the objective to protect the internal parts of the Implant, when the surgeon decides to install a temporary prosthesis after the osseointegration phase. It has a thread in its apex, connective with the respective Implant.

PRODUCT COMPOSITION:

All the Dental Implants and Cover Screws are manufactured with Degree 4 Pure Titanium in compliance with rule ASTM F67.

USE INDICATION:

The DÉRIG Dental Implants are applied in oral surgery, and are recommended for patient with unit, partial or total lack of teeth. The main function of the Implant is to serve as a support abutment to receive the fixed or removable prosthesis, in order to replace the natural root of the tooth.

The DÉRIG Implants are divided into 5 groups, and each group recommended for a region of the dental arch or in function of the class of the bone, of the technique or surgical situation adequate to each patient.

NOTE:

The **DÉRIG Cylindrical** Implants have a high power of cutting and are recommended for unit and multiple prostheses in type I, II, III and IV bones. In type I bones, the use of the Thread Forming device, available in the DÉRIG system (not provided).

The **DÉRIG Tapered** Implants have a high power of compaction and are recommended for unit and multiple prostheses in type III and IV bones.

PRECAUTIONS AND WARNINGS:

The DÉRIG Implants should be installed in controlled surgical conditions. The implants should be installed through a surgery by a qualified and experienced professional.

The professional should plan the surgery through image diagnosis tools (x-rays, tomography exams, etc) and develop a surgical plan selecting the implant(s) that are most suitable to the quantity and quality of bones available in the patient. Use of implants in patients in the growing phase is not recommended.

So that the implant can be installed correctly and integrates with the bone structure, a high accuracy in construction of the bone cavity is essential, through use of calibrated and sterile instrument and the sequence of drills corresponding to the diameter of the implant in perforation, with the adequate rotation and irrigation, thus, avoiding the thermal trauma e assuring the success of product application.

The **Dental Implant** is sterilized through gamma radiation and has a unique use and its usage is not recommended if the packing has been opened or violated prior to the moment of application. New sterilization performed by the professional is not recommended. ITS'S FORBIDDEN TO REPROCESS.

The professional should submit the patient to mouth disinfection prior to application of the product, thus preventing the implant to be in contact with a non sterile substance contaminating the bone cavity.

The professional should check the information included in the product's packing, specifically the validity term and the identification of the product.

The professional should pay attention to the torque indicated in application of the product in order not to damage it. It is recommended not to exceed the maximum torque of insertion according to the table as a guarantee of implant integrity.

Descriptive table of maximum torque applied:

Trading Model:				Maximum Torque
Cylindrical	Biodent	HEX	NP	35Ncm
Cylindrical	Biodent	HEX	RP	45Ncm
Cylindrical	Biodent	HEX	UP	45Ncm
Cylindrical	Biodent	HEX	WP	45Ncm

Trading Model:				Maximum Torque
Cylindrical	Dynamic	CMI	3.0	25Ncm
Cylindrical	Dynamic	CMI	NP	35Ncm
Cylindrical	Dynamic	CMI	RP	45Ncm
Cylindrical	Dynamic	CMI	RP	45Ncm
	Kort	HEX	NP	35Ncm
	Kort	HEX	RP	45Ncm
	Kort	CMI Plus	NP	35Ncm
Tapered	Bioneck	TRI	NP	35Ncm
Tapered	Bioneck	TRI	RP	45Ncm
Tapered	Bioneck	TRI	UP	45Ncm
Tapered	Bioneck	TRI	WP	45Ncm
Tapered	Singular	CMI	NP	35Ncm
Tapered	Singular	CMI	RP	45Ncm
Tapered	Singular	CMI	RP	45Ncm
Cylindrical	Exakort	HI	NP	45Ncm
Cover Screws (all models), Maximum torque 10Ncm <i>*see accessories</i>				

We recommend that the products be organized in such a manner that the information and shelf life term indicated in the packing can be visible.

The implants are identified according to the rules in force. The identification label provided inside the product packing should be attached to the patient medical record.

The professional should inform the patient as related to the adequate oral hygiene and the need of a periodic post-surgery follow-up.

Note: Important Precautions and Warnings for the Kort Implants.

In the region of the jaw, special attention should be given to the milling, during the osteotomy and in installation of the Kort Implant, in order to prevent the alveolar nerve to be reached or pressed as this may result in partial or total, temporary or final loss of sensibility in the affected region.

In this region, a thorough planning is needed, as well as request computerized tomography for a correct proximal sizing of the alveolar nerve, maintaining an accurate and safe distance of the nerve.

As prevention, in addition to the care mentioned above, it is recommended that upon reaching the average point of the osteotomy, to stop the action and have an x-ray taken for a new evaluation and confirmation of the plan.

CONTRAINDICATIONS:

Application of implants is not indicated in patients with vascular disorder, uncontrolled diabetes or other illnesses of metabolic or systemic nature that affect healing of tissues, coagulation disorder, treatment with anticoagulants, metabolic bone illness, patient not prepared psychologically to be submitted to full oral rehabilitation, inadequate oral hygiene, insufficient spaces between the arcs, non treatable occlusion/articulation problems, insufficient bone height and/or width, active intraoral infection and in case of abuse of alcohol and smoking. The temporary contraindications include: Chemotherapy and radiotherapy as well as chronic periodontal inflammation or insufficient coverage of soft tissue.

STORAGE, TRANSPORTATION, PRESERVATION AND MANIPULATION:

The DÉRIG implants should be stored e/or transported in their original packing and should be maintained in a dry location and temperature of 25°C \pm 2 and humidity 70%. It is recommended that the products be organized in a manner so as to allow that the information and validity terms indicated in the packing are visible.

The DÉRIG implants are identified according to the rules in force. The identification label provided inside the product packing should be attached to the patient medical record.

DISPOSAL:

All the products and material used in the dental implant installation surgery can put at risk the health of those that handled them after their usage.

Prior to disposal, we recommend checking and complying with the legislation in force

SHELF LIFE AND LOT:

Manufacturing date, shelf life and lot, see packing.

STERILIZATION:

The Implant is supplies STERILE. The DÉRIG osseointegrating implant system is sterilized by gamma rays with a charge of 25 KGy, thus assuring the sterile condition of the product, except in cases of violation of packing and products with the validity term expired.

CONTENT:

The DÉRIG implants are packed, sterile, in a close bottle with cover. The cover screws are placed, sterile, in this bottle through the cover screw holder closed with a cover. The bottled is placed in a cradle, closed through crystal cover blister and a paper of surgical degree. All of these are packed in an external cartridge. They are supplied with the following content:

- 01 Implant;
- 01 Cover Screw;
- 03 Lot Identification Labels;

HANDLING INSTRUCTION:

Check if the size corresponds to your surgical plan. Check the shelf life. Open the cardboard transportation packing. Check if the internal sterile packing has its blister perfectly sealed. Reject the product if the packing is damaged or with the shelf life expired. Stick the separate labels on the patient medical record. Open the blister only at time of implant insertion.

Start perforation with the sequence of drills corresponding to the implant.

Ask your assistant to remove the capsule within the blister. Take the bottle and position it in the cradle. Open the bottle unthreading the cover. Remove the cover screw holder to expose the implant. Capture the implant with the adequate insertion wrench. Insert the implant in the cavity recently perforated. Maximum insertion speed of 20RPM. If

done in a manual manner use the maximum torque according to the table. Open carefully the cover screw holder. Capture the cover screw with the adequate wrench. Thread it to the implant.

ADVERSE EFFECTS:

If the technique used is not the adequate technique and if the patient is not submitted to the exams duly recommended, the result of the product application can be unsatisfactory, and may cause to the patient a non-necessary re-absorption of the bone structure.

The surgical process can result in adverse effects in the region where it was applied, such as: chronic ache, tissue reaction, infection, difficulty to talk, numbness in the lip region and lower jaw in case of surgery in the lower jaw and numbness in the side tissue of the nose in case of surgery in the upper jaw. The cases of numbness are temporary effect, in very rare cases, permanent numbness occurred.

REUSE

The Dérig Implant is provided for a unique usage, thus, the re-processing is prohibited and as informed in its identification label.

Reuse of the product may result in contamination of the patient, as the sterilization is valid only for the first use of the product.

ATTENTION

**Product for dental usage, to be used only by professionals qualified for technical application.
Supplied STERILE by Gamma Rays.**

Manufacturer: Dérig Ind. e Com. de Materiais Médico-Odontológicos Ltda.

Rua Lapa, 479 - CEP 06419-020

Barueri - São Paulo - Brasil

Phone/Fax 55 (11) 4161-1991

www.derig.com.br - derig@derig.com.br

Technical Responsible Professional: Edson Aparecido Meronho –
CREA-SP nº 5063423447

ANVISA Registration N°: 80165910012

European Representative: Exaktus – Material de Reabilitação Oral

Located at Rua Laborim 68, 4430-128 Vila Nova de Gaia,
Portugal

Telefone: +351 224 068 098 | Fax +351 224 068 108

E-mail: geral@exaktus.pt

SYMBOLISM



Consult instruction for use



Do not re-use




Sterilized using radiation



Batch Code



Reference code



European Representative



Manufacturer



Date of manufacture



Use by-date



Keep protect from Sun



Don't watering



Do not use if the package is damaged



Emission: 24/01/2018

This instruction of usage was edited originally in the Portuguese language

Annex a: list of product codes.

Technical Name: Dental Implants (Osseointegrable)

Trade name: Implant

Code	Description
02.03.16.026	Cylindrical Biodent HEX NP 3.3 x 8.5mm
02.03.16.027	Cylindrical Biodent HEX NP 3.3 x 10mm
02.03.16.028	Cylindrical Biodent HEX NP 3.3 x 11.5mm
02.03.16.029	Cylindrical Biodent HEX NP 3.3 x 13mm
02.03.16.030	Cylindrical Biodent HEX NP 3.3 x 15mm
02.03.16.006	Cylindrical Biodent HEX RP 3.75 x 8.5mm
02.03.16.007	Cylindrical Biodent HEX RP 3.75 x 10mm
02.03.16.008	Cylindrical Biodent HEX RP 3.75 x 11.5mm
02.03.16.009	Cylindrical Biodent HEX RP 3.75 x 13mm
02.03.16.010	Cylindrical Biodent HEX RP 3.75 x 15mm
02.03.16.016	Cylindrical Biodent HEX RP 4,0 x 8,5mm
02.03.16.017	Cylindrical Biodent HEX RP 4,0 x 10mm
02.03.16.018	Cylindrical Biodent HEX RP 4,0 x 11,5mm
02.03.16.019	Cylindrical Biodent HEX RP 4,0 x 13mm
02.03.16.020	Cylindrical Biodent HEX RP 4,0 x 15mm
02.03.16.021	Cylindrical Biodent HEX RP 4,0 x 17mm
02.03.16.036	Cylindrical Biodent HEX UP 5.0 x 8.5mm
02.03.16.037	Cylindrical Biodent HEX UP 5.0 x 10mm
02.03.16.038	Cylindrical Biodent HEX UP 5.0 x 11.5mm
02.03.16.039	Cylindrical Biodent HEX UP 5.0 x 13mm
02.03.16.040	Cylindrical Biodent HEX UP 5.0 x 15mm
02.04.16.001	Tapered Bioneck TRI NP 3.5 x 8mm
02.04.16.002	Tapered Bioneck TRI NP 3.5 x 10mm
02.04.16.053	Tapered Bioneck TRI NP 3.5 x 11.5mm
02.04.16.003	Tapered Bioneck TRI NP 3.5 x 13mm
02.04.16.004	Tapered Bioneck TRI NP 3.5 x 16mm
02.04.16.059	Tapered Bioneck TRI NP 4.3 x 8mm
02.04.16.060	Tapered Bioneck TRI NP 4.3 x 10mm

Code	Description
02.04.16.061	Tapered Bioneck TRI NP 4.3 x 11.5mm
02.04.16.062	Tapered Bioneck TRI NP 4.3 x 13mm
02.04.16.063	Tapered Bioneck TRI NP 4.3 x 16mm
02.04.16.005	Tapered Bioneck TRI PR 4.3 x 8mm
02.04.16.006	Tapered Bioneck TRI RP 4.3 x 10mm
02.04.16.054	Tapered Bioneck TRI RP 4.3 x 11.5mm
02.04.16.007	Tapered Bioneck TRI RP 4.3 x 13mm
02.04.16.008	Tapered Bioneck TRI RP 4.3 x 16mm
02.04.16.025	Tapered Bioneck TRI UP 5.0 x 8mm
02.04.16.026	Tapered Bioneck TRI UP 5.0 x 10mm
02.04.16.055	Tapered Bioneck TRI UP 5.0 x 11.5mm
02.04.16.027	Tapered Bioneck TRI UP 5.0 x 13mm
02.04.16.028	Tapered Bioneck TRI UP 5.0 x 16mm
02.04.16.009	Tapered Bioneck TRI WP 5.0 x 8mm
02.04.16.010	Tapered Bioneck TRI WP 5.0 x 10mm
02.04.16.011	Tapered Bioneck TRI WP 5.0 x 13mm
02.04.16.012	Tapered Bioneck TRI WP 5.0 x 16mm
02.12.16.001	Cylindrical Dynamic CMI NP 3.5 x 8.5mm
02.12.16.002	Cylindrical Dynamic CMI NP 3.5 x 10mm
02.12.16.003	Cylindrical Dynamic CMI NP 3.5 x 11.5mm
02.12.16.004	Cylindrical Dynamic CMI NP 3.5 x 13mm
02.12.16.005	Cylindrical Dynamic CMI NP 3.5 x 15mm
02.12.16.006	Cylindrical Dynamic CMI RP 4.3 x 8.5mm
02.12.16.007	Cylindrical Dynamic CMI RP 4.3 x 10mm
02.12.16.008	Cylindrical Dynamic CMI RP 4.3 x 11.5mm
02.12.16.009	Cylindrical Dynamic CMI RP 4.3 x 13mm
02.12.16.010	Cylindrical Dynamic CMI RP 4.3 x 15mm
02.12.16.011	Cylindrical Dynamic CMI RP 5.0 x 8.5mm
02.12.16.012	Cylindrical Dynamic CMI RP 5.0 x 10mm
02.12.16.013	Cylindrical Dynamic CMI RP 5.0 x 11.5mm
02.12.16.014	Cylindrical Dynamic CMI RP 5.0 x 13mm
02.12.16.015	Cylindrical Dynamic CMI RP 5.0 x 15mm

Code	Description
02.12.16.101	Tapered Single CMI NP 3.5 x 8mm
02.12.16.102	Tapered Single CMI NP 3.5 x 10mm
02.12.16.103	Tapered Single CMI NP 3.5 x 11.5mm
02.12.16.104	Tapered Single CMI NP 3.5 x 13mm
02.12.16.105	Tapered Single CMI NP 3.5 x 16mm
02.12.16.106	Tapered Single CMI RP 4.3 x 8mm
02.12.16.107	Tapered Single CMI RP 4.3 x 10mm
02.12.16.108	Tapered Single CMI RP 4.3 x 11.5mm
02.12.16.109	Tapered Single CMI RP 4.3 x 13mm
02.12.16.110	Tapered Single CMI RP 4.3 x 16mm
02.12.16.111	Tapered Single CMI RP 5.0 x 8mm
02.12.16.112	Tapered Single CMI RP 5.0 x 10mm
02.12.16.113	Tapered Single CMI RP 5.0 x 11.5mm
02.12.16.114	Tapered Single CMI RP 5.0 x 13mm
02.12.16.115	Tapered Single CMI RP 5.0 x 16mm
02.03.16.049	Kort HEX NP 4.0 x 5.5mm
02.03.16.050	Kort HEX NP 4.0 x 7.0mm
02.03.16.051	Kort HEX RP 5.0 x 5.5mm
02.03.16.052	Kort HEX RP 5.0 x 7.0mm
02.17.16.001	Kort CMI Plus NP 3.5 x 5.5mm
02.17.16.002	Kort CMI Plus NP 3.5 x 7.0mm
02.17.16.003	Kort CMI Plus NP 4.3 x 5.5mm
02.17.16.004	Kort CMI Plus NP 4.3 x 7.0mm
02.17.16.005	Kort CMI Plus NP 5.0 x 5.5mm
02.17.16.006	Kort CMI Plus NP 5.0 x 7.0mm
02.15.16.001	Cylindrical Exakort HI 3.3 x 8mm
02.15.16.002	Cylindrical Exakort HI 3.3 x 10mm
02.15.16.003	Cylindrical Exakort HI 3.3 x 11.5mm
02.15.16.004	Cylindrical Exakort HI 3.3 x 13mm
02.15.16.005	Cylindrical Exakort HI 3.3 x 15mm
02.15.16.006	Cylindrical Exakort HI 3.8 x 8mm
02.15.16.007	Cylindrical Exakort HI 3.8 x 10mm

Code	Description
02.15.16.008	Cylindrical Exakort HI 3.8 x 11.5mm
02.15.16.009	Cylindrical Exakort HI 3.8 x 13mm
02.15.16.010	Cylindrical Exakort HI 3.8 x 15mm
02.15.16.011	Cylindrical Exakort HI 4,2 x 8mm
02.15.16.012	Cylindrical Exakort HI 4.2 x 10mm
02.15.16.013	Cylindrical Exakort HI 4.2 x 11.5mm
02.15.16.014	Cylindrical Exakort HI 4.2 x 13mm
02.15.16.015	Cylindrical Exakort HI 4.2 x 15mm